



BMDS

DAS-8029 XPT, IPT/T Handheld Reader

Display Only

Description

The DAS-8029 Reader is a handheld scanning device for reading ID data from BMDS XPT, IPT and IPTT implantable transponders. As a simple to use, display only reader, the DAS-8029 boasts a super bright, easy to read, OLED display with large character size. Conveniently, the last ID read will remain on the screen until a new scan or powered off.

There is an easy to follow front panel menu for selecting various option such a read beep, vibrate, and duplication filters. Auto modes including time between reads and timed shut off. Features and setup can be password protected to avoid accidental changes.

The DAS-8029 has a build in charging system using long lasting lithium rechargeable batteries. If there is ever a need to change the batteries, they are user replaceable.

The reader's control panel includes a row of sequencing LED lights that provide indication that the reader is actively scanning and attempting to capture a transponder ID. This allows the display to be momentarily shut down to minimize interference.



In an optional continuous read mode, the reader can be operated directly off it's in its AC power supply. In such use, the reader first uses the line voltage, saving the battery for portable use. In a pinch the reader will work with no batteries installed. Though the DAS-8029 is primarily a handheld system, there are special order options for stand mounting. A soft protective boot is available.





Feature

Super Bright 2 line 20 Characters per Line OLED Display

Separate Sequencing LED Indicating Active Scan providing Better Feedback to User

Hybrid Scanning Coils that Produces More Forward Read Zone.

Power Efficient Design allows All Day Usage

Control Panel

Battery Level Display Icon

Full Menu Driven Setup with Optional Password Protection

Menu Options:

- Auto Repeat Scan
- Audio Beep
- Vibrate Upon Read
- Scan Delay for Auto Mode
- Auto Mode Dup Filter
- Power Save/Timed Off
- Password for Menu
- Firmware Version

Power

12 Volt International Power Supply for Recharging

Direct Power and Battery Charging via Easy Barrel Connection

Battery: Two 3.7 Li-ion 3000mAH User Replaceable

Mechanical

Reader: 2.9" Wide by 2.0" High by 9.9" Long

Weight: 12.8 Ounces (363 Grams)

Material: Rugged ABS Polycarbonate Blend

Compatibility

BMDS XPT, IPT & IPTT Transponder Types (IPTT: No Temperature, Only ID)

CALL US: 800 526 BMDS

Bio Medic Data Systems, Inc. 1 Silas Road Seaford, DE 19973 Inquiry www.BMDS.com